

Think Twice

BETTER YOU KNOW, BEFORE YOU GO

Single-Use Breathalyzer

A cost-effective way to test for alcohol and promote workplace safety. Simply press both ends together to activate the device. After blowing into the tube and following instructions, the white crystals will change different shades of red depending on the user's breath-alcohol-content.

4-Level BAC Result

- Individually packaged
- Small, portable size
- No calibration needed
- 0.00% - 0.08%+ BAC



- ✓ FDA Approved
- ✓ Immediate Results
- ✓ Accurate & Reliable
- ✓ Millions Sold



(866) 657-5658
DUiprevention.org